



Near Patient Sharps Disposal™ & HTM 07-01

How does it work?



Point 1: Segregate Sharps & Non-Sharps Waste
Using Near Patient Sharps Disposal™ (NPSD™) at the 'point of use' limits the amount of non-sharps waste being disposed of into a sharps container, helping to reduce the cost of disposal.

Point 2: Segregate at the 'Point of Use'
Selecting the correctly coloured container and tray for the procedure means that sharps waste can be safely segregated at the 'point of use', helping to prevent needlestick injuries by eliminating unnecessary transportation of unsecured sharps.



NON-SHARPS WASTE
**Dispose of into a
WASTE SACK**
In accordance with your
waste policy.



**Dispose of into an NPSD™
container at the 'Point of Use'**
In accordance with your
waste policy.

Sharpsafe®
designed to protect

**FRONTIER
MEDICAL**