

# Near Patient Sharps Disposal™ & HTM 07-01

## How does it work?



### Point 1: Selection & Preparation

Select the correctly coloured container and tray for the procedure.

### Point 2: 'Point of Use'

Take the container and tray to the point of use. Using an NPSD™ system helps to prevent needlestick injuries by eliminating unnecessary transportation of unsecured sharps.

### Point 3: Clinical Waste Disposal

Dispose of sharps waste into the sharps container and non sharps waste into clinical waste sack in accordance with waste policy correct segregation helps to reduce the cost of disposal.

### Point 4&5: Temporary/Final Closure & Preparation

Ensure container is placed in temporary closure position in preparation for further procedures. (Engage final closure and dispose of if full) Clean tray according to local infection control policy.

