

sharpsafe[®]

1.8L container now available on drug tariff

introducing the sharpsafe[®]
1.8L container designed for
the safe vertical fill of insulin
and needle pens

sharpsafe.co.uk



helping your patients dispose of sharps safely

If your patient requires safe sharps disposal at home, look no further than the Sharpsafe range.

New to the range is the 1.8L container, designed specifically for the disposal of needle pens.

Each container is able to hold approximately 20* pens vertically, maximising the capacity of the container and preventing the need for larger containers in the patients home.



CONTACT US

	YELLOW LID Medicinal Sharps Waste		PURPLE LID Cyto Sharps Waste		
	FRONTIER PRODUCT CODE	PIP CODE	FRONTIER PRODUCT CODE	PIP CODE	DT PRICE
SHARPSAFE 1.8L	40502431	403-0342	40502421	403-0359	£1.15

01495 235 833

info@sharpsafe.co.uk

For further information on the Sharpsafe 1.8L container and the rest of our range available on Drug Tariff, please visit: www.sharpsafe.co.uk



*Based on a standard needle pen width of 17mm.