



The Sharpsafe® SMS Range

The sharps management system incorporating NPSD™ 'point of use' disposal and ANTT in a range of colours in accordance with HTM 07-01: Safe Management of Healthcare Waste.

The safer system of sharps segregation



Product features

- 1 LINE and LABEL**
Aligning the fill line on the lid with the label on the container ensures the internal safety flap inside the lid operates efficiently.
- 2 Simple FOUR CLICK assembly**
An audible click confirms the lid and base are assembled easily, correctly and securely.
- 3 Translucent lids**
A unique safety feature which indicates when the contents reach the fill line, without having to open or look into the aperture.
- 4 Clear segregation SIGNPOSTING**
✓ Waste type printed on the lid to assist with easier and safer segregation of sharps waste.
- 5 Large aperture with internal safety flap**
The large aperture with keyhole features ensures safe and easy disposal of sharps.
- 6 Temporary closure**
For increased user and patient safety when the container is not in use.
- 7 Fill line**
The inkjet printed fill line is clearly visible externally and in combination with the translucent lid and large 'headroom' ensures Sharpsafe® is the safest sharps disposal system available.
- 8 In-mould label**
This unique label is moulded directly into the container providing a permanent audit trail which cannot be removed.

Sharpsafe[®] SMS



Sharps Containers

| Products | point two | point three | point four five | point six | one | two | three | four | seven | nine | thirteen | twenty four | thirty |
|--|-------------------------|----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| INFECTIOUS SHARPS FOR INCINERATION | Cat No. 4173 NHS No. | 4175 FSL359 | 4170 FSL360 | 4150 FSL336 | 4160 FSL078 | 4140 FSL426 | 4145 FSL276 | 4100 FSL030 | 4110 FSL031 | 4117 FSL290 | 4115 FSL291 | 4120 FSL308 | 4180 FSL066 |
| INFECTIOUS SHARPS FOR INCINERATION OR ALTERNATIVE TREATMENT | Cat No. NHS No. | | | 4150B FSL350 | 4160B FSL337 | 4140B FSL159 | 4145B FSL351 | 4100B FSL338 | 4110B FSL353 | 4117B FSL341 | 4115B FSL100 | 4120B FSL342 | 4180B FSL354 |
| CYTO SHARPS FOR INCINERATION | Cat No. NHS No. | | | 4150C FSL335 | 4160C FSL321 | 4140C FSL322 | 4145C FSL277 | 4100C FSL326 | 4110C FSL283 | 4117C FSL327 | 4115C FSL329 | 4120C FSL334 | 4180C FSL279 |

Leakproof Containers



The Sharpsafe[®] Leakproof range is available in 6 Litre, 10 Litre and 19 Litre containers which comply with both UN 1851 and UN 3291.

| Products | six | ten | nineteen |
|---|----------------|--------|----------|
| INFECTIOUS SHARPS FOR INCINERATION | Cat No. 4165L | 4114L | 4119L |
| CYTO SHARPS FOR INCINERATION | Cat No. 4165CL | 4114CL | 4119CL |
| LIQUID PHARMA WASTE FOR INCINERATION | Cat No. 4165DL | 4114DL | 4119DL |

Solid Pharmaceutical Waste Containers

The Sharpsafe[®] Solid Pharmaceutical Waste container has been designed for POM's and other solid pharma waste in compliance with UN 3249.

| Products | seven | thirty |
|--|---------------------------------|-----------------|
| SOLID PHARMA WASTE FOR INCINERATION | Cat No. 4110D NHS No. FSL284 | 4180D FSL278 |

Specimen Transport Bags



| Products | STB1 6 Litre | STB2 19 Litre |
|---|--------------|---------------|
| COMMUNITY SPECIMEN TRANSPORT BAG | Cat No. 5001 | 5002 |

The Community Specimen Transport Bag (STB) has been designed for community healthcare workers who need to transport diagnostic specimens whilst adhering to UN3373.

Protected Access Containers



Designed for use in public access areas such as A&E and out-patients. The Protected Access range allows easy disposal of sharps but provides added security for both user and the general public.

The Protected Access range also discourages the use of the sharps container for non sharps waste, avoiding unnecessary waste disposal costs.

| Products | nine | thirteen | twenty four | thirty |
|--|-----------------------------------|-------------------|-------------------|-------------------|
| INFECTIOUS SHARPS FOR INCINERATION | Cat No. 4117RA NHS No. FSL292 | 4115RA FSL361 | 4120RA FSL363 | 4180RA FSL285 |
| INFECTIOUS SHARPS FOR INCINERATION OR ALTERNATIVE TREATMENT | Cat No. 4117RAB NHS No. FSL037 | 4115RAB FSL362 | 4120RAB FSL364 | 4180RAB FSL355 |
| CYTO SHARPS FOR INCINERATION | Cat No. 4117RAC NHS No. FSL328 | 4115RAC FSL366 | 4120RAC FSL365 | 4180RAC FSL356 |

Specialist Containers

| Quiver [™] | XL | Clinisafe [®] |
|---------------------|------|------------------------|
| 4190 FSL080 | 4195 | 4116 FSL079 |

FSL codes apply to NHS Supply Chain references only.

NPSD[™] System

Near Patient Sharps Disposal[™]

Sharpsafe[®] SMS, Sharps Management System is a sharps containment and management system which incorporates Sharpsafe[®] Near Patient Sharps Disposal[™] (NPSD[™]) and Aseptic Non Touch Technique (ANTT[®]).

What is NPSD[™]?

Near Patient Sharps Disposal[™] (NPSD[™]) is a strategy for reducing needlestick injuries which features sharps containment products designed specifically to be taken to the patient. Sharps are disposed of safely at the 'point of use' thus eliminating many of the opportunities for needlestick injuries.

NPSD[™] + ANTT : Aseptic Non Touch Technique


Aseptic Non Touch Technique (ANTT) is a DOH and EPIC2 endorsed 'best practice' aseptic technique developed to standardise practice, improve clinical standards and further protect patients and healthcare providers. The Sharpsafe[®] system incorporates trolleys and trays designed for use with ANTT.

Tray
Cat No. 4300



Pedestal Trolley

Cat No. 4480



Either or both




Either




Sharpsafe[®] Trolley

Cat No. 4400



Trolley & Tray System

The Sharpsafe[®] Trolley can be configured with a variety of Sharpsafe containers and/or a variety of trays.



ANTT Tray configuration for designated work areas.



Frontier can also supply a dedicated phlebotomy trolley.

Clinical Waste Containers

| | | | |
|----------|--------------|----------|--------------|
| 30 Litre | Cat No. 4185 | 60 Litre | Cat No. 4186 |
|----------|--------------|----------|--------------|



Product Stands

| | |
|------------------|----------------|
| Community Stands | Desktop Stands |
|------------------|----------------|



Brackets

| | | | |
|--------------|-----------------|--------------------------|--------------|
| Hook Bracket | Trolley Bracket | Vertical Trolley Bracket | Wall Bracket |
|--------------|-----------------|--------------------------|--------------|

